



Acousto-Optic Frequency shifters

AO optic devices frequency shift the output first order diffracted light beam by an amount equal to the RF driver (carrier) frequency. This shift will be positive or negative depending on the Bragg orientation. The frequency shift can be varied within the RF bandwidth of the AO device by tuning the input RF carrier frequency.



Maximum diffraction efficiency into the first order is achieved when the input light beam to the AO device is aligned at the Bragg angle and the RF power is adjusted for saturation.

Application Note



Frequency shift range

Isomet offer a wide range of AO modulators/frequency shifters. Typical centre frequencies range from 30 - 300MHz for devices operating in the visible to NIR and from 30 - 80 MHz for devices operating in the Far IR. However frequencies outside these general limits can be achieved using multiple AO devices or multiple passes in one AO device.

Small frequency shifts

For frequencies much lower than the minimum AO centre frequency two devices can be placed in series, one applying a frequency upshift and one a downshift. The final first order output frequency will be given by the difference of this up and down shift.


The lowest value shift is usually limited by the accuracy and stability of the driver frequency source. Values of less than 10 KHz are possible with care.

Large frequency shifts

For frequencies much higher than the maximum AO centre frequency then two or more devices can be placed in series, each applying a frequency upshift. The final first order output frequency will be the sum the upshifts.

Alternatively a multiple pass can be made through a single AO device using (say) a retro reflector. The suitability of this method will depend on the beam size relative to the maximum useable crystal aperture. The diagram below shows a basic double downshift configuration.

In the above cases, the overall efficiency will equate to the combined diffraction efficiencies. e.g. for a 85% maximum per AO modulator, the overall efficiency for a double AO system would be approximately 70%.



Application Note



Bragg Angle

The input Bragg angle is given by:

$$\theta_{\text{bragg}} = \frac{\lambda \cdot f_c}{2 \cdot v}$$

where : v = acoustic velocity
fc = carrier frequency
 λ = optical wavelength

Variable Frequency Shifts

For dynamic frequency shifting applications with fixed Bragg angle adjustment, the diffraction efficiency of the AO device will vary due to Bragg angle mis-match error. The transducer characteristics of most AO devices are designed to minimise this effect. However in practice the useable bandwidth of a single electrode AO device is limited to approximately 1/3 to 1/2 of the centre frequency.

To circumvent this difficulty and achieve greater bandwidth, the acoustic signal in the AO material can be steered and made to track the optimum Bragg conditions over a wider range of frequencies. This requires an array of electrodes on the device transducer, each with an RF signal progressively delayed. The fixed delay results in a change in phase between electrodes proportional to the drive frequency and this effects a corresponding change in the angle of the transmitted acoustic beam.

RF Drive Power

Diffraction efficiency is a function of the acoustic drive power. The optimum (saturation) power is given by:

$$P_{\text{sat}} = \frac{\lambda^2 \cdot H}{2 \cdot L \cdot M_2}$$

where: λ = optical wavelength
H = electrode height of transducer
L = electrode length
M₂ = figure of merit for interaction material

For a given device L, H and M₂ are constant and thus the value depends on the operating wavelength.

Application Note



RF Power Limitations

If the device is operated at a RF power (P) less than Psat, the reduction in efficiency can be determined from the formula:

$$\text{Relative eff' } \varepsilon = \frac{\sin^2 \frac{\pi}{2} \sqrt{P/P_{\text{sat}}}}{2}$$

This may be the case when operating at longer wavelengths with devices that have a limitation on the maximum input RF power. It is also common for the RF drive power to be tailored or programmed to flatten the amplitude of the deflected output with changing frequency.

RF Drivers

Isomet offer a comprehensive range of RF drivers. Most feature amplitude control for modulation or switching of the first order beam. The crystal controlled drivers offer the best frequency stability for fixed frequency applications, typically $\pm 0.003\%$.

For frequency tuning and frequency differential applications we recommend the iDDS-2SE dual output frequency synthesizer (plus amplifier). This offers stable and adjustable frequency differential control between to AOM's operated in series and allows precise generation of frequency shifts as low as 10KHz.

Fibre coupled Devices

Input and/or output fibre optic coupling can be provided using Fibre-optic adjuster FOA1065. The coupling efficiency of the fibre system is approx. 65-70% excluding the AO device. Four axis manipulators provide fine adjustment of the fibre optic angles. The AO mount also allows Bragg adjustment.

Typically, a fibre optic-AO frequency shifter is operated at fixed frequency. Changing frequency will necessitate re-adjustment of the AO device angle and fibre manipulators to compensate for the change of Bragg and Separation angles and thus maintain optimum coupling from the input to output fibre.